

## GST protein Isolation

**KEEP EVERYTHING ON ICE AT ALL TIMES**

### Day 0 (optional)

- Perform transformation in BL21 cells
- Plate in Lb-agar + antibiotic

### Day 1

- Start culture from transformation (Day 0) or glycerol stock
  1. Pick large single colony OR
  2. Scrape with a pipet tip over frozen glycerol stock
- Place it in 5 ml LB-medium + antibiotic
- Grow for 5-6 hours at 37 degrees in shaking incubator
- Transfer to 250 ml LB + antibiotic
- Incubate o/n 37 degrees

### Day 2

- Transfer culture 1:10 to 500 ml LB + antibiotic
- Optional - spin down rest and resuspend pellet in 50% glycerol and store at -80 degrees Celsius
- Incubate 2 hours at 37 degrees Celsius in shaking incubator
- Add 0.3 mM IPTG (Stock is 1 M - so 150 ul in 500 ml culture)
- Incubate o/n at 25 degrees Celsius in shaking incubator

### Day 3

- Spin down culture 10 min, 4000 rpm, 4 degrees Celsius (JLA 16-250 - Beckman)
- Resuspend pellet in 10 ml **Lysis buffer + Triton**
- Split suspension in two and sonicate (in 50 ml tubes in beaker with ice+bit water)
  1. Time 2 minutes amplitude 50
  2. Pulls on 7.5 seconds / pulls off 20 seconds
- (if sound changes its ok but always just finish cycle)
- Spin down 10 minutes, 12000 rpm 4 degrees Celsius (JA 30-50 - Beckman)
- Transfer supernatant to 15 ml tube and keep on ice
- For each protein prepare 500 ul of *Glutathion sepharose 4B* beads for fast flow (wash 3 times with **Lysis buffer + triton**)
- Add 500 ul of beads to each supernatant and incubate for at least 1 hr at 4 degrees rotator
- Spin down, resuspend beads in 1 ml **lysis buffer + triton** and transfer to 1,5 ml ep
- Wash 2x with **lysis buffer + triton** (30 seconds +- 3000 rpm)
- Wash 3x with **lysis buffer - triton** (30 seconds +- 3000 rpm)
- Resuspend beads in 1 ml **lysis buffer WITHOUT triton** and perform protein assay (Pierce)
- Store beads at 4 degrees Celsius

#### Pierce BCA protein assay (Kit NAME)

1. Mix 50 parts of reagens A with 1 part of reagents B
2. Add 10 ul of standard (Floris), blank, or samples in 96 well flatbottom plate

3. Add 200 ul of mixed reagents (A+B), mix 30 seconds on plate shaker/table
4. Cover plate and incubate 30 minutes at 37 degrees Celsius
5. Read absorbance at 562 nm (Magellan program -> Create/edit method (pick one) -> start measurement -> choose method -> when done view result -> transport to excel)
6. Possibly also load on SDS-PAGE and check Ponceau staining for protein presence and size

#### **LYSISBUFFER**

**+ TRITON**

**50 mM Tris, pH  
7,4**

**50 mM NaCl**

**5 mM MgCl<sub>2</sub>**

**1 mM DTT**

**1% Triton X-100**

#### **LYSISBUFFER**

**-TRITON**

**50 mM Tris, pH 7,4**

**50 mM NaCl**

**5 mM MgCl<sub>2</sub>**

**1 mM DTT**